


AL-CAN Aluminium Cylinder is the Right Choice for all your home care **oxygen Therapy, emergency medical services** (EMS), and Medical institutional Breathe easy

AMBULANCE



MOUNTAIN



ASTHMA



HOSPITAL



POLLUTION



INDUSTRIAL



SPORTS



OFFICE
(Fatigue -stress)



Available in various sizes:

In Refillable Pack - 100L to 1500L and
In Disposable pack : 10 L.

Re-Filling your Oxygen Cylinder:

The Oxygen Cylinder can be refilled with oxygen from any Medical oxygen gas manufacturer / Distributor in your area, as you will be provided the authenticate certificate issued from CCOE, Nagpur.

For Medical Gases

Water Cap.		Gas Capacity		Outside Dia		Length		Tare Weight	
Liters	C.ft.	C.mtr	Liters	Inch	mm	Inch	mm	Lbs	Kgs
1.8	9.5	0.270	270	4.48	114	12.4	315	4.0	1.8
2.2	11.6	0.330	330	4.48	114	14.0	355	4.6	2.2
3.1	16.4	0.465	465	4.48	114	18.5	470	6.8	3.1
4.5	23.8	0.675	675	4.48	114	24.8	630	9.9	4.5

Distributed by:

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313

Mob: 91 93441 66778

This seamless compressed gas cylinder is manufactured from high strength Aluminium alloy as per BS EN 1975:2000 / IS 15660:2006. The cylinder are duly inspected and Certified by TUV, India and approved by Chief Controller of explosives(CCOE), Nagpur.

Hitwin
HEALTHCARE

Life is Better with  *Pure*



- Salient Features
- Light Weight & Portable
- Corrosion Resistant
- Non-Magnetic
- Easy to handle
- Refillable
- Suitable for Hospitals,
- Nursing Homes & Home Care.