

OPTIMO COMPRESSOR NEBULIZER

MODEL: NEB 102

Hitwin
HEALTHCARE

Respiratory • Sleep Therapy • Cardiology



Highlights

- ~ Compressor Nebulizer
- ~ Durable
- ~ Reliable
- ~ Efficient Drug Deposition
- ~ Higher Dosing Accuracy
- ~ Portable

DEVICE CHARACTERISTICS

Operating Humidity Range Max(%rh)	-	95
Operating Humidity Range Min(%rh)	-	10
Operating Temperature Range Max(Degc)	-	40 Degree C
Operating Temperature Range Min(Degc)	-	10 Degree C
Storage Humidity Range Max(%rh)	-	95
Storage Humidity Range Min(%rh)	-	10
Storage Temperature Range Max(Degc)	-	70 Degree C
Storage Temperature Range Min(Degc)	-	-20 Degree C

DIMENSIONS

Length (mm)	-	300 mm
Width (mm)	-	190 mm
Height (mm)	-	95 mm
Weight (Gram)	-	1.8 g

MEDICAL SPECIFICATIONS

Average Nebulization Rate (ML/min)	-	0.4 ML/min remains amount: 0.15ml
Medication Capacity (ml)	-	6 ml
Mmad (Micro m)	-	2.44 Micro m
Particle Size (Micro m)	-	0.5 to 6 Micro m

POWER REQUIREMENTS

Power Consumption	-	below 65 W
Power	-	230 or 120V AC, 50/60 Hz
Power Source	-	Electricity

Hitwin
HEALTHCARE

Hitwin Healthcare Pvt. Ltd.
7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.
Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313
Mob: 91 93441 66778

