




Now with more flexibility

OptiLife nasal mask gives you and your patients more flexibility than ever.

With updated features and more options, it's simpler to fit more patients with a mask that's not only right for their lives ... it's a great fit for your business.

OptiLife patients can choose between a pillows cushion or a CradleCushion that literally cradles the nose to provide excellent comfort and seal. And because OptiLife is packaged with multiple cushion styles and sizes, first-time fit is straightforward, inventory management is enhanced, and callbacks may be less likely for homecare providers.

Patients will appreciate how easy OptiLife is to wear, clean, and maintain – not to mention that they have a clear field of vision and can sleep on their sides. OptiLife can be slipped into place with just one hand thanks to an advanced headgear design and an integrated chin support band. Adjustments are simple, too. Longer, multi-directional mask tubing provides increased flexibility and a secure fit, regardless of sleep position.

Hitwin
HEALTHCARE

Respiratory • Sleep Therapy • Cardiology

PHILIPS
RESPIRONICS
sense and simplicity

Easy-to-adjust headgear has no buckles – and keeps the forehead and vision area open.

Two headgear slots give options to improve patient comfort and seal.

Chin support band provides mask stability and a secure fit.

Longer, 18"/45.7 cm multi-directional tubing moves freely. The mask stays put for an optimal seal.

Simple attachment hub provides a secure cushion attachment and allows the cushion to rotate for an optimal fit.

Exhalation port provides quiet and diffused airflow.

Visual indicators provide orientation for proper cushion assembly.

Generous landing area and deep arc of the optional **CradleCushion** cradle the nose for excellent comfort and seal.



Groove accommodates side-to-side movement.

Ordering information

OptiLife mask	Pillows cushions sizes	CradleCushions sizes	Part number	U.S. HCPCS code
Mask with headgear and pillows cushions FitPack	P, S, M, L		1036800	A7034, A7035
Mask with headgear FitPack/DuoPack	P, S, M	S, M	1036801	A7034, A7035, A7032
Mask with headgear FitPack/DuoPack	M, L	M, LN, L	1036802	A7034, A7035, A7032
Mask only FitPack	P, S, M, L		1036804	A7034
Mask sub-assembly with headgear			1036818	A7035*
Mask sub-assembly			1036819	*

* Requires addition of mask cushion for HCPCS A7034

Replacement parts

Item	Size	Part number	U.S. HCPCS code
Pillows cushion, 1 per pack	P	1036838	A7033
Pillows cushion, 1 per pack	S	1036839	A7033
Pillows cushion, 1 per pack	M	1036840	A7033
Pillows cushion, 1 per pack	L	1036841	A7033
CradleCushion, 1 per pack	S	1036846	A7032
CradleCushion, 1 per pack	M	1036847	A7032
CradleCushion, 1 per pack	LN	1036848	A7032
CradleCushion, 1 per pack	L	1036849	A7032

Replacement parts

Item	Part number	U.S. HCPCS code
Headgear including chin support band	1036850	A7035
Accessory pouch	1036851	N/A
Tube management clip, 2 per pack	1039967	N/A

Respironics, CradleCushion, DuoPack, FitPack, and OptiLife are trademarks of Respironics, Inc. and its affiliates.

Please visit www.philips.com/respironics



© 2010 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Geyer JJ 11/05/10 MCI 4103751 PN 1039633

Philips Healthcare is part of Royal Philips Electronics

Philips Respironics International Headquarters
+33 1 47 28 30 82
Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9666 4444
Philips Respironics United Kingdom
+44 800 1300 845

Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668

Distributor

Hitwin Healthcare Pvt. Ltd.
7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.
Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313
Mob: 91 93441 66778