



Karma VIP 515 Tilt-in-space Wheelchair

**Redefining the concept
of comfort**



WE
ARE THEIR
WIFE LIFE
MATTERS...
Hitwin

***Folding
Tilt-in-space
Wheelchair***

Hitwin
HEALTHCARE

Respiratory • Sleep Therapy • Cardiology

Hitwin Healthcare Pvt. Ltd.

7/10, Jaganathan Nagar 2nd Main Road,
Arumbakkam, Chennai, India - 600106.

Ph: 91 44 4684 5555 (20 Lines), 91 44 4233 3313

Mob: 91 93441 66778

www.hitwinhealthcare.com



USE TILT-IN-SPACE TO IMPROVE POSITIONING

The VIP-515 is a revolutionary lightweight folding tilt-in-space wheelchair, with a tilt range from 0 to 35 degrees. With the same quality and design features of the well-known Karma 8000 series, it offers a cost effective tilt-in-space option. It is ideal for people in aged care facilities or within the community. With ergonomically designed seating that incorporates some cushioning, it can be used without additional seating products. Its contoured headrest also offers additional support. With a lightweight, folding frame it is easily transportable.

FEATURES

- Tilt-in-space from 0 to 35 degrees
- Foldable and compact
- Height adjustable and flip-back armrest
- Swingaway legrests that can rotate both inwards and outwards to allow for easy transfers
- Adjustable length legrests
- Ergonomic headrest
- Antitippers
- Heel loops on each legrest




SPECIFICATIONS



Seat Width
16", 18"



Seat Depth
17"



Overall Width
25", 27"



Overall Length
43"



Overall Height
53"



Seat to Floor Height
20"



Back Height
18"



Rear Wheel
14"



Front Castor
7"




Armrest Height
9"-11"



Total Weight
18.9 kg



User Weight Limit
115 kg



Folded Width
14"



Tilt Range
0 - 35°